



“EZ-Click” PTCA Adapter




“EZ-Click” PTCA Adapter

Main features:

- Utilizes a unique engagement arm allowing free and easy opening in three stages
- Rapid and efficient insertion and removal of any angioplasty devices
- Compatible for all known and new generation of angioplasty devices
- 3 operative stages
 - Closed
 - Partially open
 - Open
- Compatible up to 9F guiding catheters
- Minimizes back bleeding during procedure

Advantages:

- **Fast and easy transfer between open and close position**
- **Visible Device status (sealed/used)**
- **Ergonomic hand fitted design**



Device Operation

Unsealed position:




Lever holds the product at the open position

Semi-Open position:

Provides Convenience at device insertion and removal stage

Sealed position:

Press Lever - the product returns to closed position



MediSense Asia Sdn. Bhd. (714530-V)
No. 21, Bandar Menjalara, Kepong,
52200 Kuala Lumpur,
Malaysia.

Tel: +60362733068

Fax: +60362720093

Email: sales@medisenseasia.com

URL: [http:// www.medisenseasia.com](http://www.medisenseasia.com)